

How to Use a Unodose Applicator



1 Remove the protective cover with a light, counter-clockwise twist and pull while holding the barrel.



2 Turn the dosing knob, located at the bottom of the barrel, clockwise for the number of prescribed clicks. You will hear and feel a click.



Note: Leave the knob in a "lock" position, and not between click positions.



3 Confirm that the cream has exited the holes of the applicator cap.



4 Apply the cream to the skin absorption site. **Rub** with the applicator cap until the cream is evenly spread. **Replace** the protective cover with a light push and clockwise turn. You will hear and feel a click.



5 Check the scale at the blue line for the number of clicks remaining. **Refill** your prescription when the blue line reaches the Reorder label.